
2012284 Fentanyl , Urine Screen with Reflex to Quantitation

FENT RFX U

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives in ARUP Standard Transport Tubes. (Min: 1 mL)

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: Specimens exposed to repeated freeze/thaw cycles. Samples collected in tubes with additives or preservatives.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month

Reference Interval:

Effective May 17, 2021

Screen cutoff Concentration: 1 ng/mL