
2012297 Tramadol, Urine Screen with Reflex to Quantitation

TRAM RFX U

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Pharmaceutical preparation. Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: **1 month**; Frozen: **1 month**

Reference Interval:

Effective May 17, 2021

Screen Cutoff concentration **100** ng/mL