

HOTLINE: Effective May 17, 2021

2012300

Zolpidem Urine with Reflex to Quantitation

ZOLP RFX U

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine with no additives or preservatives to an ARUP Standard Transport Tube. (Min: 2 mL)

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: Samples collected in tube with additives or preservatives.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month