

HOTLINE: Effective May 18, 2020

2012572	Immunofixation Electrophoresis, Serum	IFE Q GEL
----------------	--	------------------

Performed: Sun-Sat
Reported: 1-5 days

Specimen Required: Collect: Serum separator tube.
Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Refrigerated.
Unacceptable Conditions: Plasma.
Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: 1 month

HOTLINE NOTE: There is a clinically significant charting name change associated with this test.
 Change the charting name for component 0050272, Immunofix Electrophoresis Gel from Immunofix Electrophoresis Gel to **Immunofix Electrophoresis Serum**.
 Change the charting name for component 2012601, EER Immunofix Electrophoresis Gel from EER Immunofix Electrophoresis Gel to **EER Immunofix Electrophoresis Serum**.