

Quarterly HOTLINE: Effective November 5, 2018

2013008

Periprosthetic Joint Infection (PJI) Detection (Synovasure)

SYNOVA PJI

Specimen Required: Collect: Synovial fluid in Plain Red.

Specimen Preparation: Transport 1 mL synovial fluid in the original collection tube. (Min: 0.5 mL) **Separate specimens must be submitted when multiple tests are ordered.**

Storage/Transport Temperature: Refrigerated.

Remarks: Specimen Source is required. (eg., left knee, right shoulder, etc.)

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable