

**TEST CHANGE**

Mitotane, Serum or Plasma

2013014, MITOT SP

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Plain ~~Red, Lavender~~red, lavender (K2EDTA), or ~~Pink~~pink (K2EDTA).

**Specimen Preparation:** Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP ~~Standard Transport Tube~~standard transport tube. (Min: 0.3 mL)  
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

**Transport Temperature:** ~~Refrigerated~~Frozen. Also acceptable: Room temperature or ~~frozen~~refrigerated.

**Unacceptable Conditions:** ~~Separator tubes~~Polymer gel separation tube (SST or PST)

**Remarks:**

**Stability:** Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 28 days

**Methodology:** Quantitative Gas Chromatography

**Note:**

**CPT Codes:** 80375 (~~Alt code: G0480~~)

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

~~By~~Refer to report