
2013034

Alpha Subunit, Free, Pituitary Glycoprotein Hormones (PGH)

A SUB PGH

Specimen Required: Collect: Serum separator tube (SST) or plain red.

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.25 mL) Freeze immediately.

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Frozen.

Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 24 hours; Frozen: 6 months