

HOTLINE: Effective May 18, 2020

2013115

Interleukin 17, Serum

IL17

Specimen Required: Collect: Serum Separator Tube (SST), or Plain Red.

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.4 mL)

Storage/Transport Temperature: **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered. Ship in an ARUP Standard Transport Tube.**

Unacceptable Conditions: Contaminated or heat-inactivated specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 30 minutes; Refrigerated: Unacceptable; Frozen: 1 year

Reference Interval:

Effective May 18, 2020

1.4 pg/mL or less

Note: Cytokine levels may demonstrate diurnal variation. For longitudinal comparison, it is recommended that cytokine levels be determined at the same time of day.

HOTLINE NOTE: There is a clinically significant charting name change associated with this test.

Change the charting name for component 2013113, Interleukin 17 from Interleukin 17 to Interleukin 17, Serum.

There is a numeric map change associated with this test.

Change the numeric map for component 2013113, Interleukin 17, Serum from XXXXXXXX to XXXXXX.X.