

TEST CHANGE

Cw Antigen Typing, Patient
2013504, CW AG

Specimen Requirements:

Patient Preparation:

Collect: Lavender (K2EDTA) or **p**ink (K2EDTA).

Specimen Preparation: Do not freeze. Transport 7 mL whole blood. (Min: 0.5 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Separator tubes.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable

Methodology: **Qualitative Agglutination**~~Hemagglutination (HA)~~

Note:

CPT Codes: 86905

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to report

Inserted Cells