

HOTLINE: Effective November 18, 2019

2013506

Sd(a) Antigen Typing, Patient

SDA AG

Specimen Required: Collect: Lavender (K₂EDTA) or Pink (K₂EDTA).

Specimen Preparation: **Do not freeze.** Transport 7 mL whole blood. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Separator tubes.

Stability (collection to initiation of testing): Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: Unacceptable