

Quarterly HOTLINE: Effective November 12, 2018

2013612

Infliximab or Biosimilar Activity with Reflex to Antibody

IFX DL R

Specimen Required: Patient Prep: Collect specimens before infliximab or biosimilar treatment.

Collect: Serum Separator Tube (SST).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Contaminated, hemolyzed, icteric, or lipemic specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 4 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)

Reference Interval:

Components	Reference Interval
Infliximab Activity w/Rflx to Ab	Not Detected
Infliximab Neutralizing Antibody Titer	Not Detected