

TEST CHANGE

Toxoplasma gondii Antibody, IgA by ELISA, Serum

2013890, TOXOG IGA

Specimen Requirements:

Patient Preparation:

Collect: Plain ~~r~~Red or serum separator tube~~Serum Separator Tube~~ (SST).

Specimen Preparation: Transfer 3 mL serum to an ARUP standard transport tube~~Standard Transport Tube~~. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions: Grossly hemolyzed, icteric, lipemic, and bacterially contaminated specimens.

Remarks:

Stability: Ambient: Undefined; Refrigerated: 1 week; Frozen: Indefinitely

Methodology: Semi-Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)

Performed: Varies

Reported: 5-12~~3-9~~ days

Note:

CPT Codes: 86777

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By Report