

TEST CHANGE

Prostate-Specific Kallikrein, 4Kscore

2014059, 4KSCORE

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST).

Specimen Preparation: Transfer 4 mL serum to an ARUP standard transport tube. (Min: 3 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen.

Unacceptable Conditions: Frozen serum separator tubes (SST).

Remarks: Biopsy history, ~~digital rectal exam (DRE) results~~, and clinical indication for ordering are required. Digital Rectal Exam (DRE) result should be provided, if available.

Stability: Ambient: Unacceptable; Refrigerated: 72 hours; Frozen: 1 month

Methodology: Electrochemiluminescent Immunoassay (ECLIA)

Note: 4 Kallikrein Biomarkers: Total PSA, free PSA, percent free PSA, intact PSA, and hK2. Digital Rectal Exam ~~A digital rectal exam (DRE) result is required and submissions~~ should indicate "abnormal", "normal", nodule or "unavailable". no nodule. Test should not be ordered if DRE has been performed within the last 4 days or if biopsy history is positive. DREs performed after collection of specimen are acceptable.

CPT Codes: 81539

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report

HOTLINE NOTE: There is a component change associated with this test. One or more components have been added or removed. Refer to the Hotline Test Mix for interface build information.