
2014059

Prostate-Specific Kallikrein, 4Kscore

4KSCORE

Specimen Required: Collect: Serum Separator Tube (SST).

Specimen Preparation: Transfer 4 mL serum to an ARUP Standard Transport Tube. (Min: 3 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Frozen.

Remarks: **Biopsy history must be provided at time of order.**

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 72 hours; Frozen: 1 month.