

Quarterly HOTLINE: Effective February 19, 2019

2014183	Ibuprofen Quantitative, Serum or Plasma	IBUPRO SP
---------	---	-----------

Methodology: Quantitative High Performance Liquid Chromatography/Tandem Mass Spectrometry

Specimen Required: Collect: Plain Red. Also acceptable: Lavender (EDTA) or Pink (K₂EDTA).
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.21 mL)
Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.
Unacceptable Conditions: Separator tubes
Stability (collection to initiation of testing): Ambient: 15 days; Refrigerated: 15 days; Frozen: 15 days