

TEST CHANGE

Rabies Antibody Screen by RFFIT, Serum

2014351, RABIES AB

Specimen Requirements:

Patient Preparation:

Collect: Plain ~~r~~Red or serum separator tube~~Serum Separator Tube~~ (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP standard transport tube~~Standard Transport Tube~~. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen~~Refrigerated~~. Also acceptable: Room temperature or refrigerated~~frozen~~.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 1 month

Methodology: Rapid Fluorescent Focus Inhibition

Performed: Varies

Reported: 21-31 days

Note:

CPT Codes: 86382

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report