
2014351

Rabies Antibody Screen by RFFIT, Serum

RABIES AB

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.25 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 1 month