

Quarterly HOTLINE: Effective February 19, 2019

2014521

Acetaminophen Quantitative, Urine

ACETA U

Methodology: Quantitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 1 mL urine to an ARUP Standard Transport Tube. (Min: 0.3 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Stability (collection to initiation of testing): Ambient: 1 month; Refrigerated: 1 month; Frozen: 1 month