

HOTLINE: Effective November 16, 2020

2014677

Expanded Carrier Screen by Next Generation Sequencing with Fragile X

ECS SEQ FX

Specimen Required: Collect: Lavender (EDTA).

Specimen Preparation: Transport 4 mL whole blood. (Min: 1 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated.

Remarks: Patient History form required.

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: Unacceptable