

TEST CHANGE

Expanded Carrier Screen by Next Generation Sequencing

2014680, ECS SEQ	
Specimen Requirements:	
Patient Preparation:	
Collect:	Lavender (EDTA).
Specimen Preparation:	Transport 4 mL whole blood. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	
Remarks:	Patient history form is required. <u>Patient must be at least 18 years of age, be pregnant, or be</u> <u>actively seeking/planning a pregnancy, or be a male partner of</u> <u>the pregnant patient.</u>
Stability:	Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: Unacceptable
Methodology:	Massively Parallel Sequencing / Polymerase Chain Reaction (PCR)
Performed:	Varies
Reported:	17-21 days
Note:	
CPT Codes:	81479
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
Reference Interval:	
By report	