

TEST CHANGE

Flunitrazepam and Metabolites, Urine Screen with Reflex to Confirmation/Quantitation

3000183, FLUNI URN

Specimen Requirements:

Patient Preparation:

Collect: ~~Urine~~ ~~Random urine~~

Specimen Preparation: Transfer 3 mL urine to an ARUP ~~standard transport tube~~ ~~Standard Transport Tube~~. (Min: 1.4 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen~~ ~~Refrigerated~~. Also acceptable: ~~Refrigerated~~ ~~Frozen~~.

Unacceptable Conditions:

Remarks:

Stability: Ambient: ~~2 days~~ ~~48 hours~~; Refrigerated: 2 weeks; Frozen: 3 months

Methodology: Qualitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Note: If screen is positive, then confirmation will be added. Additional charges apply.

CPT Codes: 80307; if reflexed, add 80346 (~~Reflexed Alt Code: G0480~~)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~Refer to~~ ~~By~~ report