

TEST CHANGE

Prostaglandin D2 (PG D2), Random Urine

3000240, PROST D2U

Specimen Requirements:

Patient Preparation: ~~Aspirin~~ Discontinue aspirin, indomethacin, ~~or~~ and anti-inflammatory medications ~~should be discontinued~~, if possible, at least for 48 hours prior to collection.

Collect: Urine

Specimen Preparation: Transfer 10 mL urine to ARUP standard transport tubes and freeze immediately. (Min: 5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: ~~1 month~~ 6 months

Methodology: Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)

Performed: Varies

Reported: 7-13 days

Note:

CPT Codes: 84150

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval: