
3000258

Genetic Carrier Screen, (CF, FXS, and SMA) with Reflex to Methylation

CF FX SMA

Specimen Required: Collect: Lavender (K₂EDTA). Also acceptable: Pink (K₂EDTA) or Yellow (ACD Solution A).

Specimen Preparation: Transport 5 mL whole blood. (Min: 3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma or serum. Specimens collected in sodium heparin or lithium heparin tubes. **Frozen specimens in glass collection tubes.**

Stability (collection to initiation of testing): Ambient: **72 hours**; Refrigerated: **1 week**; Frozen: **1 month**