

HOTLINE: Effective November 18, 2019

3000400 QuantiFERON-TB Gold Plus, 1-Tube

QFT-PLUS

Specimen Required: Collect: QuantiFERON-TB Gold Plus 1-tube (ARUP Supply #54015) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787. For collection and transport instructions refer to QuantiFERON under Special

Handling at https://aruplab.com/testing/specimen/quantiferon. Specimen Preparation: Transport 6 mL whole blood. (Min: 5 mL).

Storage/Transport Temperature: Refrigerated. Must be collected and shipped directly to ARUP the same calendar day.

Remarks: Do not collect or ship on, or the day before, holidays.

Stability (collection to initiation of testing): Ambient: 3 hours; Refrigerated: 48 hours; Frozen: Unacceptable