

Quarterly HOTLINE: Effective November 12, 2018

---

<b>New Test</b>	<b>3000462</b>	<b>Immature PLT Fraction</b>	<b>IMM PLT</b>
-----------------	----------------	------------------------------	----------------

---

**Methodology:** Automated Cell Count  
**Performed:** Sun-Sat  
**Reported:** Within 24 hours

**Specimen Required:** Collect: Lavender (EDTA) or Pink (K<sub>2</sub>EDTA).  
Specimen Preparation: Transport 3 mL whole blood. (Min: 0.5 mL)  
Storage/Transport Temperature: Refrigerated.  
Unacceptable Conditions: Clotted or hemolyzed specimens.  
Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 48 hours; Frozen: Unacceptable

**Reference Interval:**

Age	Male (%)	Female (%)
0-180 days	2.3-7.1	1.6-7.1
6-23 months	1.7-4.1	1.7-4.8
2-5 years	1.4-3.9	1.3-3.9
6-11 years	1.3-5.2	1.3-5.0
12-17 years	1.9-6.4	1.7-6.7
18 years and older	1.0-11.4	1.0-11.4

**CPT Code(s):** 85055

New York DOH approval pending. Call for status update.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.