

HOTLINE: Effective February 22, 2022

---

3000480

Comprehensive Systemic Sclerosis Panel

SCL COMPRE

**Specimen Required:** Collect: Serum Separator Tube (SST).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 3 mL serum to an ARUP Standard Transport Tube. (Min: 1.5 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated or contaminated specimens

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month