

Quarterly HOTLINE: Effective November 12, 2018

New Test	3000537	Eosinophil Cationic Protein (ECP)	EC PROT
-----------------	----------------	--	----------------

Methodology: Quantitative Fluorescent Enzyme Immunoassay
Performed: Sunday
Reported: 1-8 days

Specimen Required: Collect: Serum separator tube (SST).
Specimen Preparation: Collect blood by venipuncture using Serum separator tube (SST). Gently invert the tube several times. Do not shake or vortex the tube. Release ECP by clotting for 60-120 minutes at room temperature (20-24° C). Centrifuge at 1000-1300 x g for 10 minutes at room temperature. Transfer 1ml of serum to an ARUP Standard Transport Tube (Min: 0.25 mL)
Storage/Transport Temperature: Frozen, Also acceptable; Refrigerated
Unacceptable Conditions: Plasma and hemolyzed serum cannot be used
Stability (collection to initiation of testing):
 After separation from cells; Ambient: 48 hours, Refrigerated: 5 Days, Frozen: 30 Days

Reference Interval: Less than 17.8 ug/L

Interpretive Data: Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

CPT Code(s): 83520

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.