

HOTLINE: Effective July 6, 2021

---

3000967

Beryllium Quantitative, Serum or Plasma

BERYLLI SP

**Specimen Required:** Collect: Royal Blue (no additive) or Royal Blue (EDTA).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 2 mL serum or plasma to an Acid Washed Transport Vial (ARUP supply #54350) available online through eSupply using ARUP Connect or contact ARUP Client Services at (800) 522-2787. (Min: 0.6 mL)

**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions: Separator tubes.

Stability (collection to initiation of testing): Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 2 weeks