

Quarterly HOTLINE: Effective February 19, 2019

New Test	3001255	14-3-3 Protein Tau, Total, CSF	14-3-3 TAU
-----------------	----------------	---------------------------------------	-------------------



CJD Surveillance CTR Test Request Form

Methodology: Qualitative Western Blot/Quantitative Enzyme-Linked Immunosorbent Assay/Qualitative Real-Time Quaking-Induced Conversion
Performed: Varies
Reported: 7-17 days

Specimen Required: Collect: CSF

Specimen Preparation: The first 2 mL of CSF that flows from the tap should be discarded. Transfer 5 mL CSF to ARUP Standard Transport Tubes and freeze immediately. (Min: 2 mL)

Storage/Transport Temperature: Frozen.

Remarks: Completed requisition form required.

Unacceptable Conditions: Specimens exposed to more than one freeze/thaw cycle.

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: Indefinitely

Reference Interval: By report

Note: Repeat testing should be collected no sooner than 2 weeks following last encounter.

CPT Code(s): 86317; 84182; 0035U

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.