

TEST CHANGE

Prion Markers (CJD), CSF

3001255, 14-3-3 TAU

Specimen Requirements:

Patient Preparation: Patient must be 12 years of age or older.

Collect: CSF

Specimen Preparation: The first 2 mL of CSF that flows from the tap should be discarded. Transfer ~~2~~5 mL CSF to ARUP standard transport tubes and freeze immediately. (Min: ~~1~~2 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN

Unacceptable Conditions:

Remarks: Completed requisition form required. Cloudy or pink specimens may result in partial results for some components.

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: Indefinitely

Methodology: ~~Qualitative Western Blot~~/Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)/Qualitative Real-Time Quaking-Induced Conversion

Performed: Varies

Reported: 7-17 days

Note: Repeat testing should be collected no sooner than 2 weeks following last encounter.

CPT Codes: 86317; ~~84182~~; 0035U

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **January 2, 2024**