

**TEST CHANGE**

Prion Markers (CJD), CSF

3001255, 14-3-3 TAU

Specimen Requirements:

Patient Preparation: Patient must be 12 years of age or older.

Collect: CSF

Specimen Preparation: The first 2 mL of CSF that flows from the tap should be discarded. Transfer 2 mL CSF to ARUP standard transport tubes or other polypropylene tubes, taking care to avoid blood contamination from the tap. **Freeze at -20 within 20 minutes of collection, and freeze immediately.** (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN

Unacceptable Conditions:

Remarks: Completed requisition form required. Cloudy or pink specimens may result in partial results for some components.

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: Indefinitely

Methodology: Quantitative Enzyme-Linked Immunosorbent Assay (ELISA) ~~/~~  
~~/~~Qualitative Real-Time Quaking-Induced Conversion

Performed: Varies

Reported: ~~12-22~~~~7-17~~ days

Note: Repeat testing should be collected no sooner than 2 weeks following last encounter.

CPT Codes: 86317; 0035U

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report



*A nonprofit enterprise of the University of Utah  
and its Department of Pathology*

Effective Date: **December 2, 2024**