



Quarterly HOTLINE: Effective February 19, 2019

New Test	3001284	Herpesvirus 6 Antibody, IgM by IFA, Serum	HHV6 AB M
-----------------	----------------	--	------------------

Methodology: Semi-Quantitative Indirect Fluorescent Antibody

Performed: Varies

Reported: 3-7 days

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).

Specimen Preparation: Transfer 0.5 mL serum to an ARUP Standard Transport Tube. (Min: 0.1 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: 1 month

CPT Code(s): 86790

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.