

HOTLINE: Effective April 17, 2019

New Test **3001409** **Immunobullous Disease Panel, Epithelial Antibody Screening Panel** **IMBULDZPAN**

Available Now

Methodology: Indirect Fluorescent Antibody/Enzyme-Linked Immunosorbent Assay

Performed: Varies

Reported: 4-9 days

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

Reference Interval: By Report

Interpretive Data: Refer to report

Note: For specimens less than 0.5 mL, call the Immunodermatology Laboratory at (866) 266-5699.

CPT Code(s): 88346, 88350 x5, 83516 x5

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.