

3001410

Basement Membrane Zone Antibody Panel

BMZ AB PAN



Immunodermatology Required Clinical
Information Form (Serum)

Methodology: Semi-Quantitative Indirect Immunofluorescence (IIF)/Semi-Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed or lipemic specimens. Plasma

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely