
3001501

CYP2C8 and CYP2C9

2C8/2C9

Specimen Required: Collect: **Lavender** (EDTA), Pink (K₂EDTA), or Yellow (ACD Solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL) OR Transport the Saliva Collection Device.

Storage/Transport Temperature: **Refrigerated.**

Unacceptable Conditions: Plasma or serum. Specimens collected in sodium heparin or lithium heparin.

Stability (collection to initiation of testing): **Ambient:** 72 hours; **Refrigerated:** 1 week; **Frozen:** 1 month