

HOTLINE: Effective April 1, 2019

New Test **3001549** **Phosphatidylcholine Antibodies - IgG, IgM and IgA** **PHOSPHA AB**

Methodology: Semi-Quantitative Immunoassay
Performed: Varies
Reported: 3-9 days

Specimen Required: Collect: Plain Red or Serum Separator Tube (SST).
Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)
Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.
Unacceptable Conditions: Specimens transported in separator tubes.
Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 2 weeks; Frozen: 3 months

Reference Interval: By report

CPT Code(s): 83520 x3

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.