

HOTLINE: Effective December 16, 2019

New Test **3001585** **Early-Onset Alzheimer's Panel, Sequencing** **ALZ NGS**
Available Now

Methodology: Massively Parallel Sequencing
Performed: Varies
Reported: 3-6 weeks

Specimen Required: Collect: Lavender (K₂ or K₃EDTA) or Yellow (ACD Solution A or B). Peripheral blood required.
Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL)
Storage/Transport Temperature: Refrigerated.
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

Reference Interval: By report

Interpretive Data: Refer to report.
See Compliance Statement C: www.aruplab.com/CS

Note: Genes tested: *APP, PSEN1, PSEN2*

CPT Code(s): 81405; 81406

New York DOH approval pending. Call for status update.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.