

HOTLINE: Effective April 23, 2019

New Test 3001697 Carbonic Anhydrase IX by IHC CAIX IHC

Methodology: Immunohistochemistry

**Performed:** Mon-Fri **Reported:** 1-3 days

Specimen Required: Collect: Tissue.

Specimen Preparation: Formalin fix (10 percent neutral buffered formalin) and paraffin embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 5 unstained (3- to 5-micron thick sections), positively charged slides in a tissue transport kit (recommended but not required), (ARUP supply #47808) available online through eSupply using ARUP Connect or contact ARUP Client Services at (800) 522-2787. (Min: 2 slides) If sending precut slides, do not over hake

Storage/Transport Temperature: Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

<u>Unacceptable Conditions:</u> Specimens submitted with non-representative tissue type. Depleted specimens.

Stability (collection to initiation of testing): Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Note:** All stains will be handled as "Stain and Return" unless a consultation is requested. To request a consultation, submit the pathology report, all associated case materials (clinical history, blocks, slides, etc.), and the Anatomic Pathology requisition form (#32960) in place of the Immunohistochemistry Stain Form.

**CPT Code(s):** 88342

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.