

HOTLINE: Effective February 22, 2022

3001781

Extended Myositis Panel

MYOS EXT

Specimen Required: Collect: Serum Separator Tube (SST).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer three 1 mL serum aliquots to ARUP

Standard Transport Tubes: (Min: 0.5 mL/aliquot)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated or contaminated specimens.

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month