

HOTLINE: Effective December 2, 2019

---

3001798

Blastomyces Antigen Quantitative by EIA, Urine

BLASTOAG U

**Reference Interval:** Not Detected

**Interpretive Data:** The quantitative range of this assay is 1.25-200 U/mL. Antigen concentrations less than 1.25 U/mL or greater than 200 U/mL fall outside the linear range of the assay and cannot be accurately quantified.

This EIA test should be used in conjunction with other diagnostic procedures, including microbiological culture, histological examination of biopsy samples, serology and/or radiographic evidence, to aid in the diagnosis of blastomycosis.

Cross-reactivity with other endemic mycoses (Histoplasma and Coccidioides) may occur. Positive test results should be correlated with other clinical findings and relevant exposure history.

**HOTLINE NOTE:** There is a unit of measure change associated with this test.

Change the unit of measure for component 3001800, Blastomyces Antigen, Urine from ng/mL to U/ml.