
New Test **3001959** **Mitochondrial Disorders Panel (mtDNA and Nuclear Genes)** **MITO PAN**



Patient History for Molecular Genetic Testing



Informed Consent for Genetic Testing

Methodology: Next Generation Sequencing
Performed: Varies
Reported: 6-7 weeks

Specimen Required: Collect: Lavender (K₂ or K₃ EDTA). Also acceptable: Buccal Swabs.
Specimen Preparation: Transport 5 mL whole blood (Min: 2 mL) or 2 buccal swabs. (Min: 2 swabs)
Storage/Transport Temperature: Refrigerated.
Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 1 week; Frozen: Unacceptable

Reference Interval: By report

CPT Code(s): 81460; 81465; 81440

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.