
3002100	Tuberous Sclerosis Complex Panel, Sequencing and Deletion/Duplication	TSC NGS
Methodology:	Massively Parallel Sequencing	
Performed:	Varies	
Reported:	3 weeks	
Specimen Required:	<u>Collect:</u> Lavender or pink (EDTA) or yellow (ACD solution A or B). <u>Specimen Preparation:</u> Transport 3 mL whole blood. (Min: 3 mL) <u>Storage/Transport Temperature:</u> Refrigerated. <u>Unacceptable Conditions:</u> Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue. <u>Stability (collection to initiation of testing):</u> Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable	
CPT Code(s):	81405, 81406, 81407	