

HOTLINE: Effective February 18, 2020

New Test	3002106	Immunofixation, Random, Urine	U IFE
-----------------	----------------	--------------------------------------	--------------

Methodology: Qualitative Immunofixation Electrophoresis**Performed:** Sun-Sat**Reported:** 1-3 days**Specimen Required:** Collect: Random urine. Also acceptable: Urine supernatant.Specimen Preparation: Transfer one 4 mL aliquot to an ARUP Standard Transport Tube. (Min: 4 mL)Storage/Transport Temperature: Refrigerated.Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 3 weeks; Frozen: 6 months**Reference Interval:** By report**CPT Code(s):** 84156; 86335

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.