

HOTLINE: Effective February 18, 2020

New Test **3002309** **Malignancy Assessment, Pelvic Mass, Overa Plus** **OVA1 PLUS**

Methodology: Electrochemiluminescent Immunoassay (ECLIA)/Fixed-Rate-Time Nephelometry

Performed: Varies

Reported: 4-8 days

Specimen Required: Collect: Serum Separator Tube (SST).

Specimen Preparation: Transfer 2.2 mL serum to an ARUP Standard Transport Tube. (Min: 1.1 mL)

Storage/Transport Temperature: Frozen. Also acceptable: Refrigerated.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 8 days; Frozen: 9 weeks

Reference Interval: By Report

Note: Biomarkers: CA-125 II, Apolipoprotein A1 (Apo A-1), Beta-2 Microglobulin (B2M), Transferrin, and Prealbumin.

CPT Code(s): 81503

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.