

HOTLINE: Effective May 17, 2021

3002463

Connective Tissue Disease First Line Panel with Reflex

CTD PAN

Reference Interval:

Test Number	Components	Reference Interval			
0050215	Double-Stranded DNA (dsDNA) Antibody, IgG by ELISA with Reflex to dsDNA Antibody, IgG by IFA	Effective May 17, 2021			
		Test Number	Components	Reference Interval	
			Double-Stranded DNA (dsDNA) Antibody, IgG by ELISA	24 IUs or less	Negative
				25 to 30 IUs	Borderline Positive
				30 to 60 IUs	Low Positive
				60 to 200 IUs	Positive
201 IUs or greater	Strong Positive				
2002693	Double-Stranded DNA (dsDNA) Antibody, IgG by IFA (using <i>Crithidia luciliae</i>)	Less than 1:10			
2002693	Double-Stranded DNA (dsDNA) Antibody, IgG by IFA (using <i>Crithidia luciliae</i>)	Less than 1:10			
0050470	Smith/RNP (ENA) Antibody, IgG				
		29 AU/mL or less		Negative	
		30-40 AU/mL		Equivocal	
		41 AU/mL or greater		Positive	
0050085	Smith (ENA) Antibody, IgG				
		29 AU/mL or less		Negative	
		30-40 AU/mL		Equivocal	
		41 AU/mL or greater		Positive	
2012074	SSA 52 and 60 (Ro) (ENA) Antibodies, IgG				
		Test Number	Components	Reference Interval	
			SSA-52 (Ro52) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive	
			SSA-60 (Ro60) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive	
0050692	SSB (La) (ENA) Antibody, IgG				
		29 AU/mL or less		Negative	
		30-40 AU/mL		Equivocal	
		41 AU/mL or greater		Positive	
0099592	Jo-1 Antibody, IgG				
		29 AU/mL or less		Negative	
		30-40 AU/mL		Equivocal	
		41 AU/mL or greater		Positive	
0050599	Scleroderma (Scl-70) (ENA) Antibody, IgG				
		29 AU/mL or less		Negative	
		30-40 AU/mL		Equivocal	
		41 AU/mL or greater		Positive	

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 0050215, dsDNA Ab, IgG w/ Reflex to IFA Titer from no numeric map to XXXX.

There is a result type change associated with this test.

Change the result type for component 0050215, dsDNA Ab, IgG w/ Reflex to IFA Titer from Alpha to Numeric.

There is a unit of measure change associated with this test.

Change the unit of measure for component 0050215, dsDNA Ab, IgG w/ Reflex to IFA Titer from no unit of measure to IU.