

HOTLINE: Effective August 17, 2020

---

**3002479 Autoimmune Liver Disease Reflexive Panel**

**LIVER PAN**

**Note:** If F-Actin, IgG by ELISA is 20 Units or greater, then Smooth Muscle Antibody (SMA), IgG by IFA titer will be added. Additional charges apply. ANA identified by indirect fluorescence assay (IFA) using HEp-2 substrate and IgG-specific conjugate at a screening dilution of 1:80. Positive nuclear patterns reported include homogeneous, speckled, centromere, nucleolar, or nuclear dots. Positive cytoplasmic patterns reported include reticular/AMA, discrete/GW body-like, polar/golgi-like, rods and rings, or cytoplasmic speckled patterns. All positive results are reported with endpoint titers at no additional charge.

**HOTLINE NOTE:** There is a reflexive pattern change associated with this test.  
Add component 3002774, Cytoplasm Pattern to reflexive orderable 3000478