

**TEST CHANGE**

**Human Immunodeficiency Virus Type 1 (HIV-1) GenoSure Archive**

3002503, HIV GSARCH

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Lavender (~~K2 or K3~~ EDTA).

**Specimen Preparation:** Freeze immediately. Transport 4 mL whole blood in the original collection tube. (Min: 1 mL)

**Transport Temperature:** CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

**Unacceptable Conditions:** Thawed specimens.

**Remarks:**

**Stability:** Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 2 weeks

**Methodology:** Polymerase Chain Reaction (~~PCR~~)/Sequencing

**Performed:** Varies

**Reported:** ~~22-28~~<sup>10-17</sup> days

**Note:** Procedure should be used for patients with documented HIV-1 infection and undetectable viral load or low level viremia.

**CPT Codes:** 87900; 87901; 87906

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

By Report