

HOTLINE: Effective May 18, 2020

New Test	3002514	Reducing Substances - Fecal	FEC RED
-----------------	----------------	------------------------------------	----------------

Methodology: Qualitative Colorimetry
Performed: Sun-Sat
Reported: 1-2 days

Specimen Required: Collect: Stool.

Specimen Preparation: Transfer 5 g stool to an unpreserved stool transport vial (ARUP Supply #40910) available online through eSupply using ARUP Connect[™] or contact ARUP Client Services at (800) 522-2787. (Min: 1 g)

Storage/Transport Temperature: Frozen.

Unacceptable Conditions: Diapers. Stool containing barium. Specimens in media or preservatives.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 2 week

Reference Interval: Normal

Note: Normal= negative or trace. Abnormal= 1+ through 4+

CPT Code(s): 84376

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.