

HOTLINE: Effective May 18, 2020

New Test	3002537	Digitoxin Quantitative, Serum or Plasma	DIGIT SP
-----------------	----------------	--	-----------------

Methodology: Quantitative Immunoassay**Performed:** Varies**Reported:** 5-8 days**Specimen Required:** Collect: Plain Red or Lavender (K₂EDTA or K₃EDTA).Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.4 mL). **Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.Unacceptable Conditions: Separator tubes.Stability (collection to initiation of testing): Ambient: Undetermined; Refrigerated: 1 week; Frozen: 3 months**Reference Interval:** By report**CPT Code(s):** 80375 (Alt code: G0480)

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.