

TEST CHANGE

Trofile (DNA) Co-Receptor Tropism

3002570, TROFI DNA

Specimen Requirements:

Patient Preparation:

Collect: Lavender (K2 or ~~K3EDTA~~~~K3-EDTA~~).

Specimen Preparation: ~~Freeze immediately. Transport~~ ~~Transfer~~ 4 mL whole blood ~~in the original collection tube, to an ARUP Standard Transport Tube and freeze immediately.~~ (Min: 4 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN.

Unacceptable Conditions: Thawed specimens.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 2 weeks

Methodology: Polymerase Chain Reaction (~~PCR~~) / Culture

Note: Recommended for patients with undetectable viral loads.

CPT Codes: 87999

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report